



# iM3<sup>®</sup> REVOLUTION 4DC Dental X-ray



Optional floor stand

## Another iM3 Veterinary First!

- Smallest focal spot available 0.4mm - better image quality
- DC technology reduces radiation dosage compared to AC
- Wireless remote control - activate x-ray outside of room
- Approved for use in all states of Australia
- Max arm extension 230cm from wall

Simple  
wireless  
control



The Revolution 4DC features the latest technology and the full-swivel ball-joint is revolutionary in every sense. The tube head revolves freely around the sphere, so that practically any position can be achieved, including vertical. An automatic touch-sensitive lock-on/lock-off device enables effortless repositioning of the tube head between exposures.



## CR7 VET - Small, Compact Dental Scanner

- Highest resolution veterinary dental scanner available
- Flexible image plates in 7 sizes
- Largest HD plate 5.7 x 9.4cm available, suitable for extremities
- Images can be viewed throughout the clinic
- No service contracts, iM3 technical support for life of scanner

### The ultimate dental radiology package:

Revolution 4DC X-ray (wall mount) and the CR7 with size 2 & 4 plates and software

iM3 – The Veterinary Dental Company  
[www.im3vet.com](http://www.im3vet.com)

**iM3**  
The Veterinary Dental Company